

cobas[®] CT/NG Swab Sample Collection for the cobas[®] 6800/8800 Systems

cobas [®] PCR Media Uni Swab Sample Kit	P/N: 07958030190	Anorectal	Endocervical	Oropharyngeal	Vaginal
		Woven Swab	✓	X	✓

cobas [®] PCR Media Dual Swab Sample Kit	P/N: 07958021190	Anorectal	Endocervical	Oropharyngeal	Vaginal
		Woven Swab	✓	X*	✓
Flocked Swab		X	✓	X	X

*The woven swab is used to remove excess mucus from the cervical os and surrounding mucosa prior to endocervical specimen collection with the flocked swab.

