

- NewYork-Presbyterian

The below list of specimens should <u>not</u> be sent to the laboratories via the Pneumatic Tube System and should only be hand-delivered:

- All Anatomical Pathology (surgical pathology and cytology) specimens
- All Cytogenetics/FISH specimens
- All Flow Cytometry specimens
- Ammonia level blood specimens
- BALs in Lukens traps
- Blood gases, including cord bloods
- Body fluid specimens
- CSF specimens
- Lactate specimens
- Platelet Function Assay specimens
- Platelet Aggregation Assay specimens
- ROTEM (Viscoelastic Testing) specimens from ACHWN/DHK (ROTEM specimens from ICU's OK to send via PTS)