

Laboratory Urine Collection Guide

Container	Information	Container	Information	Container	Information
	<p>Tests: Drug Screen (Toxicology), HCG, Random Chemistry Urine Tests</p> <p>Collect: BD Vacutainer UA No Additive</p>		<p>Tests: Urinalysis</p> <p>Collect: BD Vacutainer UA Preservative Tube</p> <p><i>NOTE: Please fill between the min. and max lines on the tube.</i></p>		<p>Tests: Urine Culture (Microbiology)</p> <p>Collect: BD Vacutainer Urine C&S Preservative</p> <p><i>NOTE: Please fill up to the min. line on the tube.</i></p>
Container	Information	Container	Information		
	<p>Tests: Chlamydia Trachomatis and Neisseria Gonorrhoea (CT/NG PCR)</p> <p>Collect: GeneXpert Urine Specimen Collection Kit</p> <p><i>NOTE: Please fill up to the dotted line on the tube.</i></p>		<p>Tests: Trichomonas Vaginalis TMA</p> <p>Collect: Hologic Aptima Urine Specimen Collection Kit</p> <p><i>NOTE: Please fill up to the fill area on the tube.</i></p>		
	<p>Tests: Random urine tests to be sent to our reference lab.</p> <p>Collect: Sterile Container</p>		<p>Tests: 24 Hour Urine Tests</p> <p>Collect: 24 Hour Urine Container</p> <p><i>NOTE: Record the start/end date/time and volume.</i></p>		

Notes:

- This does not encompass all tubes/containers and tests that are applicable.
- Full Test requirements can be viewed on our test menu website: www.testmenu.com/nypwestchester

For any additional questions, please call the laboratory at 914-787-3276