

## Laboratory Services

### Random Urine Collection

Your doctor has ordered tests that require a random urine specimen. Below are instructions for collecting a urine specimen. It is important to follow these instructions. Any changes may result in inaccurate test results. If you have any questions, be sure to ask your provider for assistance.

#### Collection of Random Urine:

Read the instructions carefully before beginning collection and follow each of the steps to ensure you collect the correct specimen for the test.

- Use the sterile screw-capped urine container provided to you for collection.
- Complete the information requested on the container label.
  - **Label must include Last Name, First Name, Date of Birth, and Date and Time of collection.**

Smith, Johnathan 1/2/1959 8/1/2024 0700
---

---

Patient label must be complete and legible for the sample to be accepted by the laboratory.

Label must be attached to the urine container, NOT the lid.

---

- Begin urinating into the toilet (not the cup).
- Midway through urination, place the container into the stream of urine and fill the container about one half to two thirds full.
- You may finish voiding into the toilet until the bladder is empty.
- Replace the cap and tighten firmly and wipe excess urine from the outside of the container.
- Wash your hands thoroughly with soap and water.
- Transport the specimen to the laboratory within 2 hours of collection.
- If this is not possible, refrigerate and transport as soon as possible within 24 hours.
  - Urine refrigerated for more than 24 hours cannot be used for culture and will be rejected by the laboratory.