











# UCD DOPLM CHEMISTRY ALTERNATE SPECIMEN TUBES

TEST		Approved tube types for laboratory testing
<ul style="list-style-type: none"> <li>BMP</li> <li>CMP</li> <li>LIPID</li> <li>HFP</li> <li>Renal Panel</li> <li>B-HCG</li> <li>Mg, Phosphorus</li> <li>Amylase</li> <li>CRP, HSCR</li> </ul>	<ul style="list-style-type: none"> <li>CK</li> <li>GGT</li> <li>Iron</li> <li>Salicylate</li> <li>Ethanol</li> <li>Osmolality</li> <li>Lipase</li> <li>Albumin</li> </ul> <p>(includes all subcomponents of panel)</p>	 <p>PST      SST      RED TOP</p>
<ul style="list-style-type: none"> <li>TROPONIN</li> <li>PROCALCITONIN</li> <li>MYOGLOBIN</li> <li>CK-MB</li> <li>BHB</li> </ul>		 <p>PST</p>
<ul style="list-style-type: none"> <li>Complement C3</li> <li>Complement C4</li> <li>Haptoglobin</li> <li>IGG, IGA, IGM,</li> <li>Prealbumin</li> </ul>	<ul style="list-style-type: none"> <li>PSA</li> <li>Estradiol</li> <li>TSH</li> <li>T4</li> <li>Digoxin</li> <li>Ferritin</li> </ul>	 <p>SST      RED TOP</p>
<ul style="list-style-type: none"> <li>Vitamin D</li> </ul>		 <p>PST      SST</p>

# UCD DOPLM CHEMISTRY ALTERNATE SPECIMEN TUBES

<ul style="list-style-type: none"> <li>• Lithium</li> <li>• Theophylline</li> </ul>	 <p style="text-align: center;"><b>RED TOP</b></p>
<ul style="list-style-type: none"> <li>• Acetaminophen</li> <li>• Carbamazepine</li> <li>• Gentamicin</li> <li>• Methotrexate</li> <li>• Phenobarbital</li> <li>• Phenytoin</li> <li>• Tobramycin</li> <li>• Valproic acid</li> <li>• Vancomycin</li> </ul>	 <p style="text-align: center;">Sodium Heparin (Dark Green) <b>RED TOP</b></p> <p>***NOTE: There are 2 Dark Green PEDI tubes. ONLY Sodium Heparin additive is approved for testing, NOT Lithium Heparin.***</p>
<ul style="list-style-type: none"> <li>• Lactate</li> <li>• Glucose tolerance test</li> </ul>	 <p style="text-align: center;">Sodium Fluoride</p>
<ul style="list-style-type: none"> <li>• BNP</li> <li>• Ammonia</li> <li>• PTH intact</li> <li>• PTH intraoperative</li> </ul>	 <p style="text-align: center;">EDTA</p>

<ul style="list-style-type: none"><li>• HCV Screen</li></ul>	 <p>SST</p>
<ul style="list-style-type: none"><li>• ED HCV (2 tubes drawn)</li></ul>	 <p>BOTH SST WITH EDTA</p>