MICROBIOLOGY SWAB / POC / COVID TEST MEDIA GUIDE



Rectal Swab:

C. DIFF SURVEILLANCE ESBL SURVEILLANCE VRE SURVEILLANCE

CARBAPENEM RESISTANCE SURVEILLANCE

Preferred:

- BLUE label BBL CultureSwab
- Lawson #: 100816



Alternative:

- RED cap BBL DUAL CultureSwab
- Lawson # 155284



Note: *BBL* swabs are **ONLY** for rectally collected Surveillance Tests. Collect one swab per test.

CANDIDA AURIS SURVEILLANCE

Preferred:

- BD ESwab WHITE-cap
- Lawson # 146106



Alternative:

- BD ESwab **BLUE**-cap
- Lawson #: 152835



VAGINAL PANEL BY PCR

- BD Molecular Swab Collection Kit
- Lawson #: 179882



MRSA SURVEILLANCE BACTERIAL CULTURE (Aerobic Only) FUNGAL CULTURE

Preferred:

- BD ESwab WHITE-cap
- Lawson #: 146106



Alternative 1:

- Fisherbrand Fisherfinest Transport Swab
- Lawson #: 156086



Alternative 2:

- BD ESwab BLUE cap
- Lawson #: 152835



BACTERIAL CULTURE (Aerobic & Anaerobic) FUNGAL CULTURE

Note: Swabs are NOT appropriate for Anaerobic culture

Preferred:

- Anaerobic Tissue Transport Medium (ATTM)
- Lawson #: 100666



Alternative:

- Anaerobic Transport Medium (ATM)
- Lawson #: 100664



Reference: Department of Pathology & Laboratory Medicine
Updated 9/5/2024

MICROBIOLOGY SWAB / POC / COVID TEST MEDIA GUIDE



COVID and Flu A/B

COVID-2019 + FLU A/B COVID-2019 RNA, QUAL

POINT OF CARE: COVID, FLU A/B and RSV

POC LIAT SARS-CoV-2 ± FLU A/B
POC LIAT COVID PCR
POC LIAT FLU A/B + RSV

Preferred:

- BD UTM
- Lawson #'s: 100983 & 153147



Alternative:

- M4RT (Note: there are various options for M4RT)
- Lawson #'s: 154003, 154004 & 154357



SARS-COV-2 ANTIGEN

- LIAISON SARS-CoV-2
- Lawson #: 164501



POINT OF CARE RAPID STREP A

Preferred:

- BD ESwab WHITE-cap
- Lawson #: 146106



Alternative:

- BD ESwab BLUE cap
- Lawson #: 152835

