

TEST CHANGE – SPECIMEN REQUIREMENTS

Notification Date: October 31st, 2024

Effective Date: November 6th, 2024

UPDATE: [QuantiFERON-TB \(Interferon Gamma Detection \(LAB6356\)\)](#) has been removed from this list. This test continues to require a special collection kit provided by the laboratory.

Explanation: We are switching our serum collection tube from a Serum Tube (red top) to a Serum Separator Tube (SST, gold top) for Immunoserology tests. This will extend the time frame for specimen stability after separation.

Test Name:





Test ID:

ANA IgG by IFA with Reflex to ENA Antibodies and dsDNA	LAB3458
Anti-Neutrophil Cytoplasmic Antibody IgG	LAB458
Anti-nuclear Antibody, IgG	LAB147
Antimitochondrial IgG Screen	LAB513
Aspergillus Galactomannan Antigen by EIA	LAB6397
Celiac Panel	LAB822
Celiac Panel with reflex EMA	LAB821
Cyclic Citrulline Peptide (CCP)	LAB851
Cytomegalovirus (CMV) IgG and IgM	LAB957
Deamidated Gliadin Peptide (DGP) IgA and IgG	LAB725
Double-stranded DNA (dsDNA) IgG	LAB648
Epstein-Barr Antibody Panel	LAB863
F-Actin (Smooth Muscle) IgG	LAB826
Hepatitis A Panel	LAB551
Hepatitis B Core IgM	LAB549
Hepatitis B e Antibody	LAB796
Hepatitis B e Antigen	LAB908
Hepatitis B Panel	LAB3070
Hepatitis B Surface Antibody	LAB472
Hepatitis B Surface Antigen	LAB471
Hepatitis B Virus Core Antibodies (Total)	LAB1242
Hepatitis C IgG	LAB868
Hepatitis C IgG with reflex PCR	LAB6233

Vanderbilt Medical Laboratories

Herpes Simplex Virus Type 1 and Type 2 IgG	LAB507
HIV P24 Antigen and HIV-1/2 Antibody with Reflex Differentiation	LAB159
Measles IgG	LAB657
Monospot	LAB482
Mumps IgG	LAB160
Myeloperoxidase (MPO) IgG	LAB3464
Parvovirus B19 IgG and IgM	LAB800
Proteinase 3 (PR3) IgG	LAB3465
Quantiferon-TB (Interferon-Gamma Detection)	LAB6356
Rickettsia rickettsii (RMSF) IgG and IgM	LAB6398
RNP IgG	LAB3466
Rubella IgG	LAB496
Scl-70 IgG	LAB771
Smith IgG	LAB2001
Specific Extractable Nuclear Antigen (ENA) Panel	LAB3457
SSA IgG	LAB344
SSB IgG	LAB345
Tissue Transglutaminase IgA and IgG	LAB723
Tissue Transglutaminase IgA and IgG with Reflex EMA	LAB721
Toxoplasma gondii IgG and IgM	LAB501
Treponemal Pallidum Antibody	LAB1197
Varicella-Zoster Virus IgG	LAB162

Vanderbilt Medical Laboratories

Current Preferred Collection Tube	New Preferred Collection Tube
<p data-bbox="305 338 651 365">BD Vacutainer® Serum Tubes (REF 367815)</p> 	<p data-bbox="976 338 1295 365">BD Vacutainer® SST Tubes (REF 367986)</p> 
<p data-bbox="172 947 574 1014">Pediatric patients: BD Microcontainer® Blood Collection Tubes, Red Top</p> 	<p data-bbox="837 947 1430 1014">Pediatric patients: BD Microcontainer® Blood Collection Tubes, Gold Top</p> 

Recommendations:

Serum specimens for Immunoserology tests should be collected in BD Vacutainer SST Tubes or BD Microcontainer Gold Top Blood Collection Tubes.

Questions:

Please get in touch with Vanderbilt Medical Laboratories Customer Service at 615-875-5227 (5-LABS) or 800-551-5227 or visit our website: [Home | Vanderbilt Medical Laboratories \(vumc.org\)](http://Home | Vanderbilt Medical Laboratories (vumc.org))

