

DATE: \_\_\_\_\_  
 Provider Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 \_\_\_\_\_  
 Office Phone: \_\_\_\_\_  
 Office contact: \_\_\_\_\_

**For Lab Use:**  
 Order Filled By: \_\_\_\_\_ Date: \_\_\_\_\_  
 Comments: \_\_\_\_\_

**FOR COURIER PICK UP, CALL:** **805-739-3333 ext.4117 (24hrs/7 days a week)**  
**Fax supply order toll free to 844-200-0103**  
*\*Please allow 2-3 days for supply orders to be filled*  
*\*Supply orders are not considered STAT*

Qty.	Unit	Item Requested
	Each	Clinical Lab Requisition
	Each	Cytology/ Pathology Lab Requisition
	Each	Custom Lab Requisition
	Pads 100	Lab location tear pad Sp/Eng
	Each	Specimen collection guide w/photos
<b>*Required when requesting forms/include current with any changes noted</b>		
<b>Tubes</b>		
	Flat 100	Light Blue Top Sodium Citrate Tube 1.8mL
	Flat 100	Light Blue Top Sodium Citrate Tube 2.7mL
	Flat 100	EDTA Lavender Top Tube 2.0mL
	Flat 100	EDTA Lavender Top Tube 4.0mL
	Flat 100	Plain Red Tube 4.0mL
	Flat 100	SST Gold Top Tube 3.5mL
	Flat 100	Pink Top Blood Bank Tube 6mL
	Flat 100	Lithium Heparin Dark Green Top Tube 4mL
	Flat 100	Light Green Lithium Heparin PST Gel Tubes 3.0mL
	Flat 100	Gray Top Sodium Fluoride Potassium Oxalate Tube 4mL
	Set	Blood Culture Bottles Adult
	Set	Blood Culture Bottles Pediatric
<b>Pap Supplies</b>		
	Flat/25	Thin Prep Vials with Brush/Spatula Combo
	Flat/ 25 .	Thin Prep Vials with Cervical Broom
<b>Pathology</b>		
	Each	20mL formalin containers
	Each	40mL formalin containers
<b>Miscellaneous</b>		
	Pack 100	Biohazard bags regular
	Each	Biohazard bags STAT
	1	Lab Lockbox for afterhours pick ups

Qty.	Unit	Item Requested
	Each	<b>ESwab™- White Cap Regular Tip;</b> R/O MRSA aerobic, anaerobic, gram stain, genital culture Vag/Rectal-Group Beta Strep (GBS) : NO FUNGAL
	Each	<b>ESwab™- Blue Cap Flexible Mini Tip;</b> aerobic, anaerobic and gram stain genital/urethral cultures, nasopharyngeal culture <b>NO FUNGAL</b>
	Each	<b>BBL Blue Cap Double Swab</b> anaerobic/aerobic & fungal cultures (collect both swabs and send to lab)
	Each	<b>3mL Regular Tip Universal Viral Transport Media</b> flocced swab (red cap swab) Herpes Culture HSV I and II, Varicella Zoster.
	Each	<b>Aptima® Mutlitest Swab Female</b> 1. Bacterial vaginosis, candida vaginitis, trichomoniasis 2. Bacterial vaginosis, candida vaginitis, trichomoniasis, including CT/GC panel
	Each	<b>Aptima® Urine Collection</b> Male/Female Chlamydia/Gonorrhea Female Chlamydia/Gonorrhea /Trichomonas
	Each	<b>Aptima® Unisex Swab Collection</b> Chlamydia/Gonorrhea /Trichomonas (Male/Female)
	Each	<b>Dual Red Swab (Strep A)</b> Collect both swabs and send to lab
	Each	<b>3mL Flexible Mini-Tip Flocced Swab (red cap)</b> w/ Universal Viral Transport Media Influenza A/B PCR , RSV, respiratory panels, Pertussis SARS-CoV-2
	Each	<b>Ova &amp; Parasite Stool Collection</b> AlcorFix (Green top)
	Each	<b>Cary Blair for Stool Culture, green lid</b>
	Each	<b>3 oz. Sterile Specimen Container</b>
	Box 20	<b>Fit-Check (IFOBT) Patient Take Home Kit</b> with prepaid mailer & instructions (include lab order)
	Each	<b>Fit-Chek (IFOBT) Specimen Collection Tube</b> (medical provider collection -digital rectal exam)