




UCD DOPLM Tube Options for Transfusion Services

Test	Approved tube types for laboratory testing
<ul style="list-style-type: none">● Type and Screen● Crossmatch● Phenotype● A1 Lectin● Antibody Titer	<div data-bbox="771 583 852 976">A purple-capped BD Vacutainer EDTA tube with a white label. The label includes 'BD Vacutainer®', 'K2 EDTA 7.2mg', 'REF 367861', and '4.0 mL'.</div> <div data-bbox="771 1031 852 1062">EDTA</div> <div data-bbox="1042 583 1123 976">A red-capped BD Vacutainer tube with a white label. The label includes 'BD Vacutainer®', 'REF 36780', and '5.0 mL'.</div> <div data-bbox="1063 1031 1166 1062">RED TOP</div>
<ul style="list-style-type: none">● All Ortho Vision Testing (Automated Instrument Testing)● DAT● Fetal Maternal Hemorrhage (FMH)● Eluate	<div data-bbox="755 1096 836 1491">A purple-capped BD Vacutainer EDTA tube with a white label. The label includes 'BD Vacutainer®', 'K2 EDTA 7.2mg', 'REF 367861', and '4.0 mL'.</div> <div data-bbox="755 1514 836 1541">EDTA</div>